




Separation of PNGase-Released and Labeled N-Glycans by HILIC Using HALO® Glycan Column

Application Note 121-GL

Digestion of N-linked proteoglycans using PNGase F releases oligosaccharides, which can be reacted with an amine via Schiff base formation. The Schiff's base derivatives (imines) can be easily reduced to form stable amine derivatives for analysis.



Many amines have been applied for labeling glycans (Harvey, 2011, J. Chromatogr. B, 879, 1196-1225). In this application brief, procainamide was chosen because of reported improvements in ESI-MS detection (Klapoetke, et. al., 2010, J. Pharm. Biomed. Anal., 53, 315-324).

Typical Labeling Conditions:

- 1) Glycan in water (up to 10% volume)
- 2) 90+% volume of:
 - 0.4 M procainamide
 - 1M sodium cyanoborohydride in 30% glacial acetic acid/70% DMSO

12-16 hr reaction at 37°C
 SEC cleanup on Sephadex G-10 minicolumn
 Absorbance Detection @300 nm or Fluorescence with Ex 330/Em 380 nm

TEST CONDITIONS:

Column: HALO 90 Å Glycan, 2.7 μ m, 2.1 x 150 mm

Part Number: 92922-705

Mobile Phase:

A: 50 mM Ammonium formate, pH 4.45

B: Acetonitrile

Gradient: 80% B to 55% B in 25 min

Flow Rate: 0.6 mL/min

Pressure: 190 bar

Temperature: 60 °C

Detection: UV 300 nm

Injection Volume: 3.0 μ L

Sample Solvent: 70/30 ACN/water

Response Time: 0.5 sec


Data Rate: 3.3 Hz

Flow Cell: 2.5 μ L semi-micro


LC System: Shimadzu Nexera

- PEAK IDENTITIES:**
1. PAm-GlcNAc₂Man₅
 2. PAm-GlcNAc₂Man₆
 3. PAm-GlcNAc₂Man₇
 4. PAm-GlcNAc₂Man₈
 5. PAm-GlcNAc₂Man₉

STRUCTURE:



Procainamide (PAm)



A fast separation of PNGase-released and procainamide-labeled N-Glycans from Ribonuclease B is accomplished with a HALO 90 Å Glycan column.

