


CANNABIS



Isocratic Separation of 18 Cannabinoids

222-CN



PEAK IDENTITIES

- 1. Cannabidivarinic acid (CBDVA)
- 2. Cannabidivarin (CBDV)
- 3. Cannabidiolic acid (CBDA)
- 4. Cannabigerolic acid (CBGA)
- 5. Cannabigerol (CBG)
- 6. Cannabidiol (CBD)
- 7. Tetrahydrocannabivarin (THCV)
- 8. Tetrahydrocannabivarinic acid (THCVA)
- 9. Cannabinolic acid (CBNA)

TEST CONDITIONS:

Column: HALO 90 Å C18, 2.7 µm, 4.6 x 150mm

Part Number: 92814-702

Mobile Phase:

A: 20 mM Ammonium Formate, pH 2.9

B: Acetonitrile Isocratic: 76% B Flow Rate: 1.5 mL/min Pressure: 231 bar

- 10. Cannabinol (CBN)
- 11. Exo-tetrahydrocannabinol (EXO-THC)
- 12. delta 9- Tetrahydrocannabinol (D9-THC)
- 13. delta 8- Tetrahydrocannabinol (D8-THC)
- 14. Tetrahydrocannabinolic acid A (THCA-A)
- 15. Cannabichromenic acid (CBCA)
- 16. Cannabicycol (CBL)
- 17. Cannabichromene (CBC)
- 18. Cannabicyclolic acid (CBLA)

Temperature: 35 °C

Detection: UV 228 nm
Injection Volume: 4.0 μL

Sample Solvent: Methanol
Response Time: 0.025 sec

Flow Cell: 1.0 µL

System: Shimadzu Nexera X2

A HALO® C18 column is used to separate a mixture of eighteen cannabinoids, showing fast results and high resolution within critical pairs. Cannabinoids are a class of chemical compounds primarily found in the marijuana plant. Many of these compounds have been found to provide medicinal benefits such as reduction in pain and inflammation.



