# Improve HPLC Separations by Carefully Matching Column Pore-Size to Solute-Size

Richard A. Henry, Consultant 983 Greenbriar Drive State College, PA 16801 <u>HyperLC@comcast.net</u>

Stephanie A. Schuster, Sr. Chemist Advanced Materials Technology, Inc. Wilmington, DE 19810 <u>sschuster@advanced-materials-tech.com</u>



**MCF 2017 RAH** 

#### **Surface Area and Pore Structure Control Retention**



- Although the greater surface area of small-pore silica can provide more retention and resolution for small molecules, larger column pores are needed to match larger molecules.
- Full surface area of a 100Å silica is not accessible to molecules such as BSA which is only ca.
  1/2 of the average pore diameter.
- For BSA, a better choice would be ca. 400-500Å where the analyte/pore radius ratio is almost 1/10.

Angiotensin II: MW 1.0 kDa, Stokes radius < 10Å BSA: MW 66 kDa, Stokes radius 35.5Å (71Å diameter)

# Column Selection and Method Development for HPLC of Large Molecules

- High speed separations of large molecules with any particle or HPLC mode depends strongly upon selecting a particle porediameter that is "adequate" for the largest analytes in the sample.
- Giddings introduced the idea of pore-confinement, where solutes lose freedom to diffuse and interact with the surface. Both equilibrium constants (entropy) and kinetics can be impacted (1,2).
- Unger cited evidence that pore-openings should be 10 times larger than analyte diameter in order to achieve unrestricted diffusion and high performance defined as speed and resolution (3-5).
- Kirkland, et al. recently described a 1000Å solid-core particle that shows the best results yet for large protein and other analytes (6).
- Matching pores to analytes may be the first and most important step in large molecule method development as convenient experiments become readily available.

#### **Pore Diameter- A New HPLC Resolution Variable**

#### Column variables (controlled by suppliers):

- -Substrate (usually silica)
- -Bonded phase (substrate *should* have minor impact)
- Pore-size (adequately large pores needed to maximize retention, efficiency and selectivity for large molecules)

Mobile phase variables (controlled by users):

- -Solvent type: ACN, MeOH and Water for RPLC or HILIC
- -Solvent strength (% organic)
- -pH (controls ionic state of substrate and phase)
- -Additive type and concentration
- -Temperature

\* Selectivity variables originally proposed by John Dolan, LC Resources, MCF 2009.

\*\* Suppliers must be willing and able to provide full details on column variables.

### **Fundamentals of HPLC Retention & Separation**

Retention by the stationary phase:

$$V_{R1} = V_M + K_{Sorp} V_S = V_0 + V_P + K_{Sorp} V_S$$
 (1)

Separation by size exclusion:

$$V_{R2} = V_0 + K_{Size} V_P$$
 (2)

 $K_{Sorp}$  can vary from 0 to infinity, and  $K_{Size}$  can vary from 0 for a large molecule to 1 for a small molecule.

These equations are usually employed separately when either the retention mode or the exclusion mode dominates.

### **Combining Equations When Modes Overlap**

When both large and small molecules are present they cannot all access pores equally, therefore both equations must be considered to interpret retention.

$$V_{\rm R} = V_0 + K_{\rm Size} \left( V_{\rm P} + K_{\rm Sorp} \, V_{\rm S} \right) \tag{3}$$

- Equation (3) more completely describes retention and separation and simply applies the pore-size constant to everything that goes on within the pores, where most of the stationary phase interaction occurs.
- For the simple case of small molecules (K<sub>Size</sub> = 1), retention equation (3) reduces to equation (1). Maximum retention is possible for small molecules where K<sub>Size</sub> = 1 and K<sub>Sorp</sub> > 0.
- Lightest retention at the excluded volume is observed for large molecules (larger than particle pores) where K<sub>Size</sub> = 0.
- Separating and analyzing large and small molecules together in a single column may be the most challenging HPLC experiment.

#### Learning About Size From SEC Calibration Curves



column: Aqueous SEC columns, 30 cm x 7.8 mm I.D., 5 μm particles mobile phase: 150 mM phosphate buffer, pH 7.0 flow rate: 1.0 mL/min det.: UV at 214 nm inj.: 10 μL

- Note that molecules 8, 9 with MW ≤1,350 can move freely within particle pores.
- Molecules 6, 7 with MW ≥10,000 are excluded from about half of the pore volume of a 150Å column.
- Molecules 1-5 require wide-pore columns for retention modes (RP, etc.)

\* Data provided by Supelco

#### **Small Molecule Retention Behavior**



\* Data provided by Supelco

**MCF 2017 RAH** 

#### **Large Molecule Retention Behavior**



#### MCF 2017 RAH

\* Data provided by Supelco

#### **Retention Decreases as Large Molecule MW Increases**

Note that 300Å C18 shows better overall retention for small and large molecules when they are present in the same sample. Widepore columns should be evaluated for this application.



Molecular weights for various proteins

**MCF 2017 RAH** 

\* Data provided by Supelco

#### **Biomolecules Often Exceed Particle Pore Diameters**

- Biomolecules are important to living organisms. They can range in size from amino acids or small lipids to large proteins such as antibodies or polynucleotides such as DNA or RNA.
- New HPLC/UHPLC product designs must rely on knowledge of how pore size and shape affects separation of large and small molecules



### **Size Comparison of Biomolecules to Pore Openings**



10 nm = 100 Å

While small pore particles may be adequate for drugs and small peptides, larger 300-1000 Å pores are needed for proteins and large biomolecules.

### Effect of Radius Ratio (a/r) on SEC Retention



Cylindrical pore model

Pore radius (**r**) varies within a range for porous particles; an average is usually reported.

Analyte radius (**a**) increases with MW; a specific ratio (**a/r)** exists for each sample component (7).

As solutes become larger, they are excluded from pores to create the calibration curve, but lose access to surface area for retention.



Conical pore model



#### **Theoretical Curve for Pure Size Exclusion HPLC**



$$V_{M} = V_{0} + V_{P} \qquad (1)$$
$$V_{R} = V_{0} + K_{Size}V_{P} \qquad (2)$$

K<sub>Size</sub> varies from zero for total pore exclusion to one for total pore occupation. Separation by a pure SEC mechanism correlates with solute hydrodynamic radius and pore geometry (8).

The subject of this study is to identify the point on the curve where solute size begins to significantly impact column performance.

**MCF 2017 RAH** 

#### **Calibration Curve for Size Exclusion Chromatography**

#### **Chromatographic Conditions**

CO	lumns: S	Sepax Zenix SEC-150, 30 cm x 4.6 mm I.D., 3 µm					
mobile	e phase: 0	.2 M potassiu	um phos	sphate, pH	7.0		
flo	w rate: 0	0.25 mL/min					
pre	essure: 6	66 bar					
column temp.: 25 °C							
de	tector: U	UV 215, 280 nm					
injection:		0.5 μL				Solute/pore	Solute/pore
samples:		listed below				radius 0 1-1 0	$\sim$ radius < 0.1
					( )		
			Size		6 -	1	
	Analyte	MW (Da)	(Å)	K <sub>size</sub>		2,3	y = -1.4279x + 8.2024
1	thyroglobul	in 667000	200	0.00	5 -	5 - 4 5	$R^2 = 0.9641$
2	SigmaMa	150000		0.06		6 7	
3	lgG	150000	100	0.05	≥ 4 -	• • • • • • • • • • • • • • • • • • • •	8
4	BSA	66400		0.16	Mode overlap	9 10	
5	ovalbumir	45000		0.24	<u>6</u> 3-		•
6	myoglobir	17000	40	0.33		region shows lower column performance	11
7	ribonuclease	e A 13700		0.38	2 -		
8	bovine insu	lin 5700		0.60			
9	neurotensi	n 1700		0.66	1 -		

0.85

0.74

1.00

•••• • 12 1  $V_{\rm M}$ V<sub>0</sub> 0 -1.50 2.00 2.50 3.00 3.50 4.00 4.50

**Elution Volume** 

Data provided by AMT

10

11

12

neurotensin

vitamin B12

angiotensin II

uracil

1700

1350

1000

112

20

#### **Impact of Size Exclusion on Band Broadening**

Sources of solute zone spreading in packed-bed liquid chromatography have been widely studied. Schure (6) recently described another type of band broadening by estimating lower diffusion rates inside crowded pores ( $D_p$ ) relative to free diffusion in mobile phase outside pores ( $D_m$ ) and plotting it against radius ratio (molecule/average pore).



#### **Two Important Porous Particle Designs**

Totally Porous Particle (TPP)





**Superficially Porous Particle (SPP)** 



Shell with 90 Å pores



#### **Examples of Halo Silica Surface-Layer Porosity**



- Unique speed and efficiency of superficially porous particles arise from the ability to use larger particles with shorter diffusion paths at lower pressure.
- Large molecules need larger pores and diffuse 10-100 times slower than small ones, benefitting more from the short diffusion paths of SPP columns.
- SPP advantages for interactive modes (RP, HILIC, NP) increase for large molecules and may extend to SEC mode; pores up to 1000Å are, available.

#### MCF 2017 RAH

### **Pore-Sizes of Halo Silica for Different Molecules**



MCF 2017 RAH

### **Pore-Crowding Lowers Efficiency for High MW Solute**

Column: 4.6 x 100 mm Halo Peptide ES-C18; T= 60 °C Sample: Bovine insulin, MW= 5733; Inj vol = 10 μL Mobile Phase: 31.5% ACN/68.5% Water/0.1% TFA, Flow Rate: 0.5 mL/min



High efficiency is the key to high peak capacity; only the larger pore particle will have high peak capacity for bovine insulin.

From a previous slide, we see that insulin has a diameter of about 50Å.

MCF 2017 RAH

# Halo 90 Å: Pores Too Small for Large Peptides

Columns: 100 mm x 4.6 mm HALO C18 (90 Å pores) and 100 mm x 4.6 mm HALO ES-C18 (160 Å pores); mobile phase: A: water/0.1% trifluoroacetic acid; B: acetonitrile/0.1% trifluoroacetic acid; gradient: 25–42% B in 10 min; flow rate: 1.5 mL/min; temperature: 30 ° C; detection: 215 nm; Peak widths in minutes above each peak.



# HALO 160 Å: Pores Too Small for Large Proteins



**MCF 2017 RAH** 

# Light Chain MAb Fragment on Halo 160Å and 400Å



Columns: 150 mm x 2.1 mm; mobile phase: A: water/0.1% difluoroacetic acid; B: acetonitrile/0.1% difluoroacetic acid; gradient: 27-37% B in 20 min; flow rate: 0.4 mL/min; temperature: 80 ° C; injection volume: 2 µL; detection: 280 nm

#### \* Data provided by AMT

#### MCF 2017 RAH

# Comparison of 160Å and 400 Å on Large Proteins

![](_page_23_Figure_1.jpeg)

MCF 2017 RAH

## Humira<sup>®</sup> (adalimumab): Too Large for 300Å Columns

![](_page_24_Figure_1.jpeg)

**MCF 2017 RAH** 

25

# SigmaMAb: Too Large for 300Å Columns

![](_page_25_Figure_1.jpeg)

**MCF 2017 RAH** 

# Small Molecules Confirm Relative Pore-Sizes and Surface Areas for 300Å and 400Å Columns

![](_page_26_Figure_1.jpeg)

### Conclusions

- Column retention and efficiency can both be lost for larger sample components when size exclusion mode significantly overlaps with retention modes and partially exclude large molecules and interfere with normal diffusion processes.
- Guidelines have been proposed that solutes should be smaller than ca. 10% of the column exclusion limit for optimum performance (radius ration of 0.1). For example, if exclusion limit is ca. 150,000 Da, largest solute size to avoid significant peak retention and efficiency loss might be estimated at 15,000 Da.
- Preliminary sample screening by SEC can save valuable time by identifying solutes that might be subject to excessive pore exclusion and require evaluation with larger pore columns.
- If data on solute sizes and average column pore size is not available to select optimum radius ratios in advance, screening larger pore columns is recommended for unknown samples.

#### References

- 1. J. C. Giddings, Unified Separation Science, Chapter 2, Entropy Effects in Porous Media, Wiley (1991).
- 2. J. C. Giddings, Dynamics of Chromatography: Principles and Theory, CRC Press, Boca Raton, FL (2002).
- 3. K. Unger, Porous Silica, Chapter 9 Size-Exclusion, Elsevier (1979).
- 4. F. Eisenbeis and S. Ehlerding, Kontakte 1/78 (1978) 22-29.
- 5. G. Barka and P. Hoffman, J. Chromatogr., 389 (1987) 273-278.
- B. M. Wagner, S. A. Schuster, B. E. Boyes, T. J. Shields, W. L. Miles, M. J. Haynes, R. E. Moran, J. J. Kirkland and M. R. Schure, J. Chromatogr. A, 1489 (2017) 75-85.
- 7. J. Klein and K. Treichel, Chromatographia, 10 (1977) 604-610.
- 8. High Performance Liquid Chromatography, J. H. Knox, ed. Edinburgh University Press (1978).

L.R. Snyder, J.J. Kirkland, J.W. Dolan, Introduction to Modern Liquid Chromatography, 3rd ed., John Wiley & Sons, Hoboken, NJ, 2010.

#### **Acknowledgements and Trademarks**

Thanks to Stephanie Schuster, Joe DeStefano and Jack Kirkland of Advanced Materials Technology for data and helpful discussions; Mark Schure also contributed helpful information.

Thanks also to William Betz, Hillel Brandes, Stacy Squillario and Roy Eksteen of Supelco Division of MilliporeSigma who also provided data.

www.advanced-materials-tech.com www.mac-mod.com

Halo and Fused-Core are registered trademarks of Advanced Materials Technology.