Advantages of Superficially Porous Particles (SPP) for Faster HPLC and UHPLC Separations

Stephanie J. Rosenberg¹, Stephanie A. Schuster¹, Conner McHale¹, Richard A. Henry² ¹Advanced Materials Technology, Inc., Wilmington, DE 19810 ²Independent Consultant, Penn State University, State College, PA 16801


🦚 advancedmaterialstechnology ted at HPLC 2019 Poster: 448

Introduction


- . Parallel efforts by James Jorgenson who developed sub-2µm, fully porous particles (FPP)¹ and Jack Kirkland who went in a new direction with superficially porous particles (SPP)² now allow much faster HPLC experiments to be carried out.
- Low, flat van Deemter plots for the latest FPP and SPP column technology show large efficiency improvements and opportunities to develop fast, new QC methods and also revalidate methods that already employ HPLC without the high pressure
- 2 um SPP column technology delivers resolution and speed advantages for optimized UHPLC and LC-MS methods.

The Unique Superficially Porous Particle (SPP)






Porous silica structure throughout



The purpose of improving efficiency is to maintain or increase resolution. Small particle sizes, such as 2 µm with shorter lengths deliver efficiency quickly, making them an ideal choice for optimized high throughput analysis. Notice that the separation above is x4 faster on the SPP particle compared to the FPP, saving not only time but solvent.



The 2 µm HALO BiohenvI is an ideal choice for high throughput analysis of drug panels, in which isobaric species separation is needed. Note the resolution between codeine and hydrocodone, (peaks 5 and 8, respectively) and morphine and hydromorphone (peaks 1 and 3, respectively).

Optimization Consideration to Achieve the Best 2 µm Performance

Optimized systems are necessary to see the most benefit from high efficiency columns, specifically SPP's of smaller diameter such as the 2 um

The following considerations are the most common and easiest to implement:

System Volume · Minimize excess tubing and install appropriately sized id tubing for the UHPLC and/ or LCMS system

Improper column dimension for LC instrument

Improper injection volume and solvent

- If using an LC detector, such as UV, choose a small flow cell volume (<1 μl ideally)
- Consider inherent LC system volumes resulting from the pumping configuration (binary vs.
- quaternary vs. isocratic), injector design (fixed loop vs. in-flow) Data Collection Considerations:
 - LC Detector Sampling Rate- 20 Hz minimum, higher preferred
 - MS Scan Speed- as fast as possible considering targeted or untargeted analysis

· Overlooked items that compromise column efficiency

Improper pore sizeImproper connections

Conclusions


· 2 μm particles deliver speed and increased efficiencies with the use of modern and optimized UHPLC and LC-MS systems

- Due to SPP performance at lower operating back pressures compared to FPP, longer column lengths can be used when high resolution is required.
- · Short SPP columns can be easily adopted for fast, high throughput LC-MS assays of complex mixtures where the MS resolving power can be exploited
- Development of new SPP columns and technologies continues with new particle geometries, pore sizes and bonded phase chemistries to further the scope of their application.
- 1. J. W. Jorgenson, et. al., LCGC N. A., Vol. 21, Number 7, July 2003, 600-610.
- 2. J. J. Kirkland, et. al., American Laboratory, 39 (February 2007), 18-21.

HALO * and Fused-Core * are registered trademarks of Advanced Materials Technology, Inc



Single process technique subject to more variability





Performance Plots for 5, 2.7, and 2 µm SPP Columns

2000

- van Deemter Equation
- H = height equivalent to theoretical plate A = eddy diffusion term
- B = longitudinal diffusion term
- C = resistance to mass transfer term μ = mobile phase linear velocity (L/t_o)
- $H = A + \frac{B}{-} + Cu$

