

LC-MS Method Development and Column Screening for Pharmaceutical and Personal Care Products (PPCPs) in the Environment

Conner McHale
Advanced Materials Technology Inc., Wilmington, DE



ASMS 2021: WP-076

Introduction

Common pharmaceutical and personal care products (PPCPs), such as over-the-counter medications, veterinary prescriptions, soaps, lotions, and even insect repellents, have become a growing concern to our environment. PPCPs can be detected at low levels in municipal wastewater, polluted ground water, and even drinking water and come from a variety of different sources. LC-MS method development and LC column screening is performed based on EPA 542 for PPCPs to determine whether phase chemistries other than the method recommended C18 will provide improved chromatographic resolution and a more accurate analysis.




Figure 1: Major pathways of PPCP release into the environment. (nih.gov)

EPA 542: Determination of PPCPs in Drinking Water by LC/ESI-MS/MS

Table 1. HPLC Conditions (Positive ion electrospray)

Time (min)	%5 mM ammonium acetate in 10% MeOH/90% reagent water*	%MeOH
0.00	90	10
0.50	90	10
0.51	50	50
8.00	25	75
8.01	0	100
10.00	0	100
15.00	90	10
24.00	90	10

Preparation of 5 mM ammonium acetate in 10% MeOH/90% reagent water: Combine 385 mg ammonium acetate and reagent water in a 1 L volumetric flask. Add 100 mL MeOH and dilute to volume.

Table 2. HPLC Conditions (Negative ion electrospray)


Time (min)	%5 mM ammonium acetate in 10% MeOH/90% reagent water*	%MeOH
0.00	90	10
0.50	90	10
0.51	40	60
8.00	0	100
11.00	0	100
15.00	90	10
25.00	90	10

Preparation of 5 mM ammonium acetate in 10% MeOH/90% reagent water: Combine 385 mg ammonium acetate and reagent water in a 1 L volumetric flask. Add 100 mL MeOH and dilute to volume.

Figure 2: EPA 542 using a HALO 90 Å C18, 2.7µm 2.1x150mm

EPA 542 recommends using two different gradients with ammonium acetate and methanol as the mobile phases. The column of choice is a C18 column. Chromatogram shows the separation of 12 PPCP compounds performed on a HALO 90 Å C18 column.

Hydrophobic Subtraction Model/ Phase Screening




- H measures the hydrophobicity of the phase.
- S* measures the resistance of the stationary phase to penetration by a solute molecule.
- A measures the hydrogen-bond acidity of the phase.
- B measures the hydrogen-bond basicity of the phase.
- C measures the interaction of the phase with ionized solute molecules and is measured at pH 2.8 and 7.0.

Table 1. HSM values

Fs	Phase	USP type	H	S*	A	B	C (pH 2.8)	C (pH 7.0)
0	HALO C18	L1	1.100	0.040	0.000	-0.050	0.050	0.040
10.04	HALO C8	L7	0.910	0.020	-0.130	0.000	-0.010	0.180
12.07	HALO AQ-C18	L1	1.000	-0.036	0.099	-0.048	0.156	0.864
17.35	HALO Phenyl-Hexyl	L11	0.780	-0.090	-0.230	0.000	0.100	0.450
17.43	HALO C30	L62	0.938	-0.046	-0.140	0.023	0.170	0.350
22.78	HALO ES-CN	L10	0.566	-0.110	-0.344	0.021	0.126	1.150
26.76	HALO Biphenyl	L11	0.708	-0.183	-0.279	0.028	0.047	0.990
52.83	HALO RP-Amide	L60	0.850	0.080	-0.380	0.190	-0.410	0.310
94.45	HALO PFP	L43	0.702	-0.073	-0.062	1.170	0.972	

*Note HSM values listed are for 2.7 µm particle size.




Column screening was performed using 9 different stationary phases (2.1x100mm) in order to increase resolution between PPCP compounds. Screening gradients were used (10-100%B in 10 min) with Water/ 0.1% formic acid and Acetonitrile/ 0.1% formic acid, 0.3 mL/min at 30°C. The HALO RP-Amide stationary phase showed best overall resolution/ peak shape.


EPA 542: Compound Table

Analyte	Drug Category	Analyte	Drug Category
1 Carbamazepine	Anticonvulsant	7 Gemfibrozil	Cholesterol Medication
2 Diazepam	Anxiolytic and Sedative	8 Naproxen	NSAID
3 Diclofenac	NSAID	9 Phenytoin	Anticonvulsant
4 Enalapril	ACE Inhibitor Medication	10 Sulfamethoxazole	Antibiotics
5 Erythromycin	Antibiotics and Gut Motility Stimulator	11 Triclosan	Antibacterial and Antifungal Agent
6 Fluoxetine	Selective Serotonin Reuptake Inhibitor	12 Trimethoprim	Antibiotics and Folate Synthesis Inhibitor

DryLab® Optimization




PPCP Analysis Using RP-Amide



DryLab® method optimization further increased peak resolution using the RP-Amide stationary phase. Methods were combined into one gradient in order to increase throughput. Gradient: 10-100% B in 24 min with Water/ 0.1% formic acid and Acetonitrile/ 0.1% formic acid, 0.3ml/min at 34°C.


Structures



RP-Amide Stability

TEST CONDITIONS:
Column: HALO 90 Å RP-Amide, 2.7µm, 2.1x100mm
Part Number: 92812-607
Mobile Phase A: Water
Mobile Phase B: ACN
Flow Rate: 1 mL/min
Temperature: 30 °C
Detection: UV, 254 nm
Back Pressure: 518 bar
Injection Volume: 0.5 µl
Sample: uracil, phenol, 1 Cl-4- nitrobenzene, naphthalene

A HALO 90 Å RP-Amide column shows excellent stability using a mixture of neutral compounds.



Conclusions

Over the next century, the combination of increasing global population size and potential droughts may result in reduced water availability, increased need for water reuse, and increasing concentrations of PPCP in wastewaters. The current wastewater treatment methods do not remove all PPCP effectively. This, coupled with the possibility that antibiotics may promote the development of antibiotic-resistant bacteria and antibiotic-resistant genes, leads to concerns about the sustainability of global water supplies. Improved chromatographic methods can help this environmental concern by showing higher resolution and efficiency compared to existing validated methods.

*HALO is a registered trademark of Advanced Materials Technology, Inc