





Improved Signal Intensity for Trastuzumab Peptides Using a HALO® 1.5 mm ID Column

BIO-325



This app note features the full MS scan of a sample of Trastuzumab that underwent trypsin digestion to produce peptide fragments of varying length that can then be used for peptide mapping. With the help of a computer program, an extracted ion chromatogram (XIC) of the sequence YADSVKGR is also featured. This XIC is to give a closer look at the benefits of a 2µm 1.5 mm ID peptide column. Samples with similar complexity require the increased efficiency from a smaller particle size and long shallow gradients. By switching from a 2.1 mm ID to a 1.5 mm ID not only can solvent usage be cut in half, but sensitivity can also be increased.

This all can be achieved by using the HALO® 1.5mm ID 2µm ES-C18 product.



**TEST CONDITIONS:**

Column: HALO 160 Å ES-C18 , 2.0 µm, 1.5 x 150 mm

Part Number: 9112X-702

Column: HALO 160 Å ES-C18 , 2.0 µm, 2.1 x 150 mm

Part Number: 91122-702

Mobile Phase A: Water, 0.1% DFA

Mobile Phase B: Acetonitrile, 0.1% DFA

| Gradient: Time | %B |
|----------------|------|
| 0.5 | 2 |
| 60.5 | 50 |
| 61.0 | 70 |
| 65.0 | 70 |
| 65.5 | 2 |
| 70.0 | Stop |

Flow Rate: 0.2 mL/min for 1.5 mm
0.4 mL/min for 2.1 mm

Pressure: 372 bar 1.5 mm
670 bar 2.1 mm

Temperature: 60 °C

Injection Volume: 1 µL

Sample: 1mg/mL Trastuzumab tryptic digest

Sample Solvent: 1.5M Guanidine HCl/0.5% Formic Acid/~50mM Tris pH: 7.8

LC System: Shimadzu Nexera X2

TUBING OPTIMIZATION:

50µm x 600mm Column to Diverter Valve

50µm x 350mm Diverter Valve to Ground

50µm x 100mm Ground to Source

MS CONDITIONS:

System: ThermoFisher Q Exactive

Spray Voltage (kV): 3.8

Capillary temperature: 320 °C

Sheath gas: 35

Aux gas: 10

RF lens: 50

