





Increased Peak Capacity of Trastuzumab Tryptic Digest on PCS C18

357



#	Tryptic Peptide	XIC	t _R (min)
1	AEDTAVYYC(Carbamidomethyl)SR	667,7877 Z=2	ES-C18: 6.41 PCS C18: 4.60
2	TPEVTC(Carbamidomethyl)VVVDVSHEDPEVK	713,6807 Z=3	ES-C18: 12.28 PCS C18: 10.11
3	TVAAPSVFIFPPSDEQLK	973,5171 Z=2	ES-C18: 17.12 PCS C18: 14.47



DIGESTION PROCEDURE:

Standard digest conditions were used for an overnight digestion of trastuzumab drug product (stock concentration 21mg/mL) at 37 °C with shaking (final concentration 1.25µg mAb/µL). The sample buffer was 50mM ammonium bicarbonate. The sample was diluted with 50mM ABC to 1.5M Guanidine prior to trypsin digestion. The next day, the digest was adjusted to 0.5% formic acid prior to LCMS analysis. 2% ACN was added to the sample prior to analysis to aid solubility. The injected sample consisted of 1.5M guanidine HCl, 2% ACN, 0.5% Formic Acid, ~50mM Ammonium Bicarbonate, 1.25µg/µL digested mAb, 0.06µg/µL trypsin.

A separation of Trastuzumab tryptic digest is performed on two 160 Å HALO columns, the ES-C18 and PCS C18 phases. On the MS system a formic acid mobile phase is used in order to maintain high ionization efficiencies. Because of the use of a low ionic mobile phase additive (formic acid) separation of the digested mAb is difficult for the standard ES-C18 phase. By using a positively charged stationary phase (PCS C18) with low ionic conditions allows for an alternative selectivity and better separation of the peptides. By measuring peak width @ 50% of 3 distinct peptides it can be seen how the effect of the PCS C18 phase can significantly help peptide separations that require low ionic mobile phases such as formic acid.



**TEST CONDITIONS:**

Column: HALO 160 Å ES-C18 , 2.7 µm, 2.1 x 150 mm

Part Number: 92122-702

Column: HALO 160 Å PCS C18 , 2.7 µm, 2.1 x 150 mm

Part Number: 92112-717

Mobile Phase A: Water + 0.1% Formic Acid

Mobile Phase B: Acetonitrile + 0.1% Formic Acid

Gradient:	Time	%B
	0.0	3
	30.0	50
	30.1	95
	33.0	95
	33.1	3
	37.0	3

Flow Rate: 0.4 mL/min.

Pressure: 465 bar

Temperature: 60 °C

Injection Volume: 1 µL

Sample: Trastuzumab Tryptic Digest (1.25 µg/µL)

Sample Solvent: Refer to Digestion Procedure

LC System: Shimadzu Nexera X2

MS CONDITIONS:

System: QExactive HF

ESI positive polarity

300-2000 *m/z*

Source voltage: 3.2kV

Sheath Gas: 40

Aux Gas: 20

Aux Gas Temp: 275 °C

Capillary Temp: 320 °C

µscans: 1

Max Injection Time: 200 msec.

S-Lens RF: 50

Tubing Optimization:

Column outlet to Diverter Valve: AMT MarvelXACT™ PEEKsil™ 50 µm ID x 350 mm

Part Number: PS7050350

Diverter Valve to Ground: AMT MarvelXACT™ PEEKsil™ 50 µm ID x 350 mm

Part Number: PS7050350

Ground to Source: AMT MarvelXACT™ PEEKsil™ 50 µm ID x 150 mm

Part Number: PS7050150

