

Drugs of Abuse and Metabolite Screening Using High Efficiency LC Columns

Stephanie Schuster, Conner McHale, Peter Pellegrinelli
Advanced Materials Technology, Inc.





advancedmaterialstechnology

Pittcon


February 28, 2024

- Background
- Superficially Porous Particles (SPPs)
- HALO® Column Chemistries used in the Screenings
- Results
- Summary



- Goal to develop an LC-MS/MS method containing fentanyl and xylazine
- Motivation:
 - Significant increase in fentanyl overdose deaths where xylazine is detected from January 2019 to June 2022 in the US
 - April 2023: White House Office of National Drug Control Policy designated fentanyl adulterated with xylazine as an emerging threat

- Horse tranquilizer that is not Schedule I restricted
 - In June 2023, Pennsylvania added it to Schedule III making it a controlled substance
- No reversal agent for humans
- Known as “tranq” or “tranq dope”
- Central nervous system depressant that can cause the following:
 - Drowsiness
 - Amnesia
 - Slow breathing, heart rate, and blood pressure to dangerously low levels
- Increases the risk of wounds that will not heal, which could lead to amputations and even death if left untreated



know what's in your drugs **TRANQ | XYLAZINE**

Xylazine is a veterinary tranquilizer that is cut in dope to give fentanyl longer legs. It's known as "anestesia de caballo" in Puerto Rico and "tranq" in Philly.

Tranq was found in over
90%
of dope samples tested in
Philly in 2021

Dope with tranq was first seen in Puerto Rico. Today, it is being found in more and more places across the US.



When tranq is mixed with another drug (like fentanyl, heroin, or a benzo), the chance of overdose increases. If someone is overdosing administer naloxone like you normally would. If the person starts breathing again but is still sedated, they don't need more naloxone. Put them in rescue position and keep an eye on them.



Tranq has been associated with severe wounds, which spread and worsen very quickly.



These wounds are seen regardless of how people use: smoking, snorting, or injecting. It's very difficult for these wounds to heal on their own so it is important to get medical attention if you have one.

What can you do if you think there is tranq in your dope?

First, try to ask around and see how the drug is making other people feel before you buy or use it. Since tranq can cause a really heavy nod, try to use somewhere that you will be safe and won't fall and hurt yourself. Finally, if you think there is tranq in your dope let others know--including someone at your local exchange program--so folks know to be careful.



carry narcan
(naloxone)



start low and
go slow




tell someone
you're using

Never Use Alone


English hotline: 800-484-3731
Spanish hotline: 800-928-5330

The Brave App

download in the app store



Created by the Substance Use and Harm Reduction (SUPHR) division
For more information on tranq in Philadelphia visit substanceusephilly.com/tranq



23 Compounds Included in Screening

Opioid Polysubstance Mix Kit	Description
6-Acetylmorphine	metabolite of heroin
Alpha-hydroxyalprazolam	metabolite of alprazolam
Alprazolam	drug of abuse
(±)-Amphetamine	drug of abuse
Benzoyllecgonine	metabolite of cocaine
Cocaethylene	metabolite formed when cocaine and ethanol are taken together
Cocaine	drug of abuse
Diazepam	drug of abuse
Diphenhydramine HCl	drug of abuse
Ecgonine HCl	metabolite of cocaine
Fentanyl	drug of abuse
Heroin	drug of abuse
Hydrocodone	drug of abuse
(-) Levamisole HCl	cocaine adulterant
(±)-Methadone	used to treat opioid addiction
(±)-Methamphetamine	drug of abuse
Morphine	drug of abuse
Naloxone	used to treat opioid overdose
Noroxycodone HCl	metabolite of oxycodone
Oxycodone	drug of abuse
Oxymorphone	metabolite of oxycodone
Tetrahydrocannabinol ((-)-delta9-THC)	drug of abuse
Xylazine	fentanyl adulterant

Opioid Polysubstance Mix Kit -
standards from Cerilliant
Xylazine - Supelco

23 Compounds Included in Screening

Opioid Polysubstance Mix Kit	Description
6-Acetylmorphine	metabolite of heroin
Alpha-hydroxyalprazolam	metabolite of alprazolam
Alprazolam	drug of abuse
(±)-Amphetamine	drug of abuse
Benzoyllecgonine	metabolite of cocaine
Cocaethylene	metabolite formed when cocaine and ethanol are taken together
Cocaine	drug of abuse
Diazepam	drug of abuse
Diphenhydramine HCl	drug of abuse
Ecgonine HCl	metabolite of cocaine
Fentanyl	drug of abuse
Heroin	drug of abuse
Hydrocodone	drug of abuse
(-) Levamisole HCl	cocaine adulterant
(±)-Methadone	used to treat opioid addiction
(±)-Methamphetamine	drug of abuse
Morphine	drug of abuse
Naloxone	used to treat opioid overdose
Noroxycodone HCl	metabolite of oxycodone
Oxycodone	drug of abuse
Oxymorphone	metabolite of oxycodone
Tetrahydrocannabinol ((-)-delta9-THC)	drug of abuse
Xylazine	fentanyl adulterant

Opioid Polysubstance Mix Kit -
standards from Cerilliant
Xylazine - Supelco

12 drugs of abuse

23 Compounds Included in Screening

Opioid Polysubstance Mix Kit	Description
6-Acetylmorphine	metabolite of heroin
Alpha-hydroxyalprazolam	metabolite of alprazolam
Alprazolam	drug of abuse
(±)-Amphetamine	drug of abuse
Benzoyllecgonine	metabolite of cocaine
Cocaethylene	metabolite formed when cocaine and ethanol are taken together
Cocaine	drug of abuse
Diazepam	drug of abuse
Diphenhydramine HCl	drug of abuse
Ecgonine HCl	metabolite of cocaine
Fentanyl	drug of abuse
Heroin	drug of abuse
Hydrocodone	drug of abuse
(-) Levamisole HCl	cocaine adulterant
(±)-Methadone	used to treat opioid addiction
(±)-Methamphetamine	drug of abuse
Morphine	drug of abuse
Naloxone	used to treat opioid overdose
Noroxycodone HCl	metabolite of oxycodone
Oxycodone	drug of abuse
Oxymorphone	metabolite of oxycodone
Tetrahydrocannabinol ((-)-delta9-THC)	drug of abuse
Xylazine	fentanyl adulterant

Opioid Polysubstance Mix Kit -
standards from Cerilliant
Xylazine - Supelco

7 metabolites

Compounds Included in Screening

Opioid Polysubstance Mix Kit	Description
6-Acetylmorphine	metabolite of heroin
Alpha-hydroxyalprazolam	metabolite of alprazolam
Alprazolam	drug of abuse
(±)-Amphetamine	drug of abuse
Benzoyllecgonine	metabolite of cocaine
Cocaethylene	metabolite formed when cocaine and ethanol are taken together
Cocaine	drug of abuse
Diazepam	drug of abuse
Diphenhydramine HCl	drug of abuse
Ecgonine HCl	metabolite of cocaine
Fentanyl	drug of abuse
Heroin	drug of abuse
Hydrocodone	drug of abuse
(-) Levamisole HCl	cocaine adulterant
(±)-Methadone	used to treat opioid addiction
(±)-Methamphetamine	drug of abuse
Morphine	drug of abuse
Naloxone	used to treat opioid overdose
Noroxycodone HCl	metabolite of oxycodone
Oxycodone	drug of abuse
Oxymorphone	metabolite of oxycodone
Tetrahydrocannabinol ((-)-delta9-THC)	drug of abuse
Xylazine	fentanyl adulterant

Opioid Polysubstance Mix Kit -
standards from Cerilliant
Xylazine - Supelco

2 adulterants

Compounds Included in Screening

Opioid Polysubstance Mix Kit	Description
6-Acetylmorphine	metabolite of heroin
Alpha-hydroxyalprazolam	metabolite of alprazolam
Alprazolam	drug of abuse
(±)-Amphetamine	drug of abuse
Benzoyllecgonine	metabolite of cocaine
Cocaethylene	metabolite formed when cocaine and ethanol are taken together
Cocaine	drug of abuse
Diazepam	drug of abuse
Diphenhydramine HCl	drug of abuse
Ecgonine HCl	metabolite of cocaine
Fentanyl	drug of abuse
Heroin	drug of abuse
Hydrocodone	drug of abuse
(-) Levamisole HCl	cocaine adulterant
(±)-Methadone	used to treat opioid addiction
(±)-Methamphetamine	drug of abuse
Morphine	drug of abuse
Naloxone	used to treat opioid overdose
Noroxycodone HCl	metabolite of oxycodone
Oxycodone	drug of abuse
Oxymorphone	metabolite of oxycodone
Tetrahydrocannabinol ((-)-delta9-THC)	drug of abuse
Xylazine	fentanyl adulterant


Opioid Polysubstance Mix Kit -
standards from Cerilliant
Xylazine - Supelco

2 treatments

Superficially Porous Particles




HALO 90 Å, 2.7 μm




HALO® Advantages

- Fast separations
- Lower back pressure compared to sub-2- μm particle columns
- Increased efficiency/resolution



HALO® Column Chemistries Selected for Screenings




HALO® Biphenyl

Features and Benefits

- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Enhanced retention and selectivity for polar molecules


HALO® Column Chemistries Selected for Screenings



HALO® Biphennyl

Features and Benefits

- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Enhanced retention and selectivity for polar molecules



HALO® RP-Amide

- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Alternate selectivity to alkyl phases

HALO®



- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Enhanced retention and selectivity for polar molecules




- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Alternate selectivity to alkyl phases

HALO® PCS C18

- Optimized for low ionic strength mobile phases
- Improved peak shape and sample loading capacity for basic compounds


HALO® Column Chemistries Selected for Screenings



HALO® Biphenyl


Features and Benefits

- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Enhanced retention and selectivity for polar molecules




HALO® RP-Amide

- Resistant to dewetting, making it 100% aqueous mobile phase compatible (can start gradients at 0% organic)
- Alternate selectivity to alkyl phases



HALO® PCS C18

- Optimized for low ionic strength mobile phases
- Improved peak shape and sample loading capacity for basic compounds



HALO® LPH-C18


- Optimized for small molecule reversed-phase HPLC separations using low pH mobile phases and elevated temperature for acidic and neutral compounds
- Unendcapped

Drugs of Abuse/Metabolites Screening with HALO® Phases



- Stationary phases (all 2.7 μm ; 90 Å in 2.1 x 100 mm) MS Conditions


- Biphenyl
- RP-Amide
- PCS C18
- LPH-C18




- Method Conditions

- Flow Rate: 0.5 mL/min
- Mobile Phase A: water/0.1% formic acid
- Mobile Phase B: methanol/0.1% formic acid
- Gradient: 1 min. hold at 3% B; 3-100% B in 9 min; 1 min hold at 100% B
- Temperature: 35 °C
- Injection: 5 μL of 50 $\mu\text{g/mL}$
- Sample: Mix of 23 drugs and metabolites
- Sample Solvent: Mobile Phase A
- Instrument: Shimadzu Nexera and LCMS 8040


- Spray Voltage (kV): 4.5
- Nebulizer Gas Flow: 2 L/min
- DL Temperature: 250 °C
- Heat Block Temperature: 400 °C
- Drying Gas Flow: 15 L/min



Screening Results: Biphenyl




HALO®




Name	Retention Time (min)
Ecgonine HCl	1.00
Morphine	3.17
Noroxycodone HCl	3.29
Oxymorphone	3.29
(±)-Amphetamine	3.42
(±)-Methamphetamine	4.03
Naloxone	4.22
6-Acetylmorphine	4.40
Oxycodone	4.53
(-) Levamisole HCl	4.61
Hydrocodone	4.65
Xylazine	5.45
Heroin	5.65
Benzoylcegonine	5.76
Cocaine	5.78
Cocaethylene	6.22
Diphenhydramine HCl	6.51
Fentanyl	6.66
(±)-Methadone	7.54
Alpha-hydroxyalprozalam	8.42
Alprazolam	8.81
Diazepam	9.18
Tetrahydrocannabinol ((-)-delta9-THC)	9.98

Screening Results: RP-Amid




HALO®




Name	Retention Time (min)
Ecgonine HCl	0.39
Morphine	0.72
(±)-Amphetamine	1.46
(-) Levamisole HCl	1.47
(±)-Methamphetamine	1.74
Naloxone	2.11
Oxycodone	2.22
Oxymorphone	2.37
Noroxycodone HCl	2.37
Hydrocodone	2.47
6-Acetylmorphine	2.78
Xylazine	3.27
Cocaine	3.59
Benzoylcegonine	3.69
Heroin	3.75
Cocaethylene	4.15
Fentanyl	4.66
Diphenhydramine HCl	4.81
(±)-Methadone	5.61
Alpha-hydroxyalprazolam	6.34
Alprazolam	6.53
Diazepam	7.10
Tetrahydrocannabinol ((-)-delta9-THC)	9.61

Screening Results: PCS C18




Name	Retention Time (min)
Ecgonine HCl	0.39
Morphine	0.60
Noroxycodone HCl	0.73
(±)-Amphetamine	1.36
(-) Levamisole HCl	1.42
(±)-Methamphetamine	1.71
Naloxone	2.15
Oxycodone	2.25
Oxymorphone	2.29
Hydrocodone	2.50
6-Acetylmorphine	2.76
Xylazine	3.31
Cocaine	3.60
Heroin	3.78
Benzoylcegonine	4.16
Cocaethylene	4.16
Fentanyl	4.75
Diphenhydramine HCl	4.81
(±)-Methadone	5.68
Alpha-hydroxyalprozalam	6.80
Alprazolam	7.01
Diazepam	7.57
Tetrahydrocannabinol ((-)-delta9-THC)	9.94

Screening Results: LPH-C18




Comparison of Retention Times



- Biphenyl was most retentive for the majority of compounds
- LPH-C18 was the second highest retentive phase
- For RP-Amide and PCS C18, retention was similar for some compounds and higher with PCS C18 for other compounds such as alprazolam and diazepam

Comparison of Tailing Factors




- Best tailing observed with Biphenyl and PCS C18


Xylazine and Fentanyl Peak Shapes: PCS C18 and LPH-C18



PCS C18




Xylazine




Fentanyl

LPH-C18




Xylazine




Fentanyl

Xylazine and Fentanyl Peak Shapes: Biphenyl and RP-Amide

Biphenyl




Xylazine




Fentanyl

RP-Amide




Xylazine




Fentanyl

Xylazine Calibration Curve




- Calibration curve for xylazine was run using the HALO® Biphenyl column
- Good linearity was observed from 10-500 ng/mL




Observed Fronting Peak Shape

- 0.1% formic acid has a pH of 2.7
- The pH is too close to the pKa of alprazolam. It should be far enough away for the alpha-hydroxyalprazolam. Both compounds show fronting.




Alpha-hydroxyalprazolam
pKa = 4.9




Alprazolam
pKa = 2.4

HALO® Biphenyl
0.1% Formic acid




Alpha-hydroxyalprazolam




Alprazolam

Improved Peak Shape with Increased pH

0.1% Formic acid/pH 2.7




Alpha-hydroxyalprazolam




Alprazolam

20 mM Ammonium Formate/ 0.1% Formic acid/pH 3.6



Alpha-hydroxyalprazolam



Alprazolam


Gradient Optimization with HALO® Biphenyl




- Faster Method
 - Removed the 1 minute hold at 3% B
 - Used a linear gradient from 0-100% B in 6 minutes
 - Kept a 1 minute hold at 100% B from 6-7 minutes
 - Total time with re-equilibration is 12 minutes
- Used the same flow rate of 0.5 mL/min
- Added 20 mM ammonium formate/0.1% formic acid to both mobile phase A and mobile phase B for improved peak shape
 - Needed to use 95/5 ACN/water for mobile phase B for solubility of ammonium formate

Faster Gradient on HALO® Biphenyl


HALO®



Impact of Reduced Flow Rate



Larger Peak Areas with 0.25 mL/min compared to 0.5 mL/min



- Increased ionization efficiency with 0.25 mL/min flow rate

- HALO® Biphenyl showed the best compromise of retention and peak shape of the 4 phases screened
- Xylazine calibration showed a linear response
- 100% aqueous compatibility enables the gradient to begin at 0% B
- HALO® columns can be run at fast flow rates for increased throughput or at reduced flow rates for increased response, depending on the requirements of the analysis

Questions

support@advanced-materials-tech.com



Photo by [Alexander Grey](#) on [Unsplash](#)



advancedmaterialtechnology



halocolumns.com



3521 Silverside Road, Suite 1-K
Quillen Building
Wilmington, DE 19810



(302) 992-8060

ISO 9001 QMS certified
All operations in Wilmington,
DE
Global distribution

