




Advantage of Inert Hardware with HALO® OLIGO C18

390



TEST CONDITIONS:

Column: HALO 120 Å OLIGO C18, 2.7 μm , 2.1 x 50 mm

Part Number: P2A62-402

Mobile Phase A: 100mM TEAA, Adjusted to pH = 8.4

Mobile Phase B: ACN

Gradient:	Time	%B
	0.0	8
	3.0	10
	3.5	20
	4.0	20
	4.1	8
	8.0	8

Flow Rate: 0.5 mL/min

Back Pressure: 146 bar

Temperature: 60 °C

Injection: 1.0 μL

Sample Solvent: 10mM Tris HCl/ 1mM EDTA pH = 8.0

Wavelength: PDA, 254 nm

Flow Cell: 1 μL

Data Rate: 40 Hz

Response Time: 0.05 sec.

LC System: Shimadzu Nexera X2

PEAK IDENTITIES

1. Oligo dT, 15 mer

Oligonucleotides are known to exhibit non-specific adsorption to stainless steel. In this comparison, the advantages of the inert column hardware over stainless steel hardware are demonstrated. The peak area is 46% larger and the tailing factor is 26% lower with the inert column hardware, which is used for HALO® OLIGO C18. Furthermore, the retention time is decreased since non-specific metal interactions are reduced when using the inert hardware.