




## LCMSMS Investigation of Tirzepatide with HALO® PCS C18

392



### TEST CONDITIONS:

Column: HALO 160 Å PCS C18, 2.7µm, 2.1 x 100mm

Part Number: 92112-717

Mobile Phase A: H<sub>2</sub>O + 0.1% Formic Acid

Mobile Phase B: Acetonitrile + 0.1% Formic Acid

Gradient:	Time	%B
	0.0	20
	2.0	20
	10.0	80
	12.0	80
	12.1	20
	15.0	20

Flow Rate: 0.3 ml/min

Pressure: 95 Bar @ 20%B

Temperature: 60 °C

Detection: PDA 280nm ; MS<sup>1</sup> and MS/MS

Injection Volume: 1 µl

Sample Solvent: 10mM Tris pH 7.8 in H<sub>2</sub>O

Data Rate: UV 40Hz

LC System: Shimadzu Nexera X2

MS System: Thermo Q-Exactive HF

### PEAK IDENTITIES

1. Tirzepatide

### MS CONDITIONS:

Ion Mode: Positive Electrospray

Sheath Gas Flow: 35

Aux Gas Flow: 10

Sweep Gas: 2

Spray Voltage: 4 kV

Capillary Temperature: 320 °C

Aux Gas Heater Temp: 275 °C

S-Lens RF Level: 60%

m/z Scan Range: 350-2000

MS1 Resolution: 120,000

DDA MS/MS Resolution: 30,000

AGC Target: 3x10<sup>6</sup>

Quad Isolation Window: m/z 2.0

Stepped NCE: 20, 30, and 40

GLP-1 agonists represent a rapidly growing market in biopharmaceuticals. These compounds represent a challenge for both chromatography and MS analysis due to their size and composition, particularly the use of lipidation to improve uptake and metabolism. This application demonstrates the capability of the HALO® PCS C18 column to maintain excellent peak profiles when weak ion pairing agents such as formic acid are required for LC/MS analysis. The HALO® PCS C18 represents an excellent choice for LCMS characterization of biologic pharmaceuticals.