



Current LC/MS Approaches for PFAS Analysis with Short and Long Chain Mixtures

Conner McHale, Barry Boyes, Joshua McBee



Conner McHale:
Technical Support Specialist
Advanced Materials Technology

cmchale@advanced-materials-tech.com
Phone: 1-302-992-8060 *1124

Presentation Outline

- Advanced Materials Technology

- Superficially Porous Particles (SPP) vs. Fully Porous Particles (FPP)
- HALO Enviroclass (PFAS Solution)

- Short and Long Chain PFAS

- Existing EPA Methods
- The Challenges
- Charge Surface RP/AEX







Image credit: NIEHS.



- Method Development

- Column Selection
- Mobile Phase Optimization
- Dry Lab Optimization



- Applications

- Well Water Samples
- Soil Samples (Collaboration w/ PFAS Solutions)



- Technical Resources/ Support

Founded in 2005 by Tim Langlois and Joe DeStefano

First company to commercially manufacture sub 3 μm superficially porous particles – *Fused-Core®*


Facility

- ISO 9001 QMS certified company
- Fully equipped state of the art laboratories
- Own process: seed reaction
- All operations handled in Wilmington, DE
 - R&D, Applications, QA/QC, Manufacturing, Sales and Marketing
- Global distribution





AMT is a company of innovators and continues to grow and deliver enabling materials to market. Our incredible team is our greatest resource.

Superficially Porous Particle (SPP)




Fully Porous Particle (FPP)




- Higher Efficiency/ Peak Capacities: SPPs provide a 30-50% increase in efficiency over FPPs of the same particle size
- Faster Analysis Times: Due to enhanced mass transfer
- Great Performance at Lower Back Pressures: SPPs allow for performance similar to smaller FPPs but at significantly lower back pressures

HALO® Enviroclass: PFAS Solution


HALO®






The PFAS Family

Panieri E, Baralic K, Djukic-Cosic D, Buha Djordjevic A, Saso L. PFAS Molecules: A Major Concern for the Human Health and the Environment. *Toxics*. 2022; 10(2):44. <https://doi.org/10.3390/toxics10020044>



PFAS Environmental Distribution and Exposure Routes

Panieri E, Baralic K, Djukic-Cosic D, Buha Djordjevic A, Saso L. PFAS Molecules: A Major Concern for the Human Health and the Environment. *Toxics*. 2022; 10(2):44. <https://doi.org/10.3390/toxics10020044>



HALO
ENVIRONMENTAL

EPA method 1633 PFAS

Chromatogram Data:


Peak #	Compound	Retention Time	Peak #	Compound	Retention Time	Peak #	Compound	Retention Time
1	PFBA	1.23135	16	ADONA	6.13967	31	NEFOSSAA	9.28502
2	PFMBA	2.46145	17	6-2FTS	6.89270	32	11C1-PF3OLHS	9.65908
3	3-3 FTCA	3.59342	18	PFDA	6.89485	33	PFDOA	9.76475
4	PFHxA	3.67722	19	PFHpA	6.92040	34	NMFOA	10.47117
5	PFBS	4.04783	20	PFNA	7.49722	35	NMFOSE	10.50321
6	PFMBA	4.13198	21	PFOS	7.72580	36	PFDoS	10.52672
7	PFESA	4.60998	22	7-3 FTCA	7.88195	37	PFHxA	10.57155
8	NFOHA	4.96555	23	PF3ONS	8.19448	38	NEFOSE	10.76262
9	4-2FTS	4.96621	24	8-2FTS	8.49818	39	NEFOA	10.76952
10	PFHxA	5.06168	25	PFDA	8.51965	40	PFHDA	11.07977
11	PFHNS	5.25443	26	PFNS	8.87143			
12	HPFO-DA	5.37952	27	NMFOSSAA	9.00733			
13	PFHDA	6.03142	28	PFOSA	9.26438			
14	PFHNS	6.11412	29	PFDA	9.28270			
15	5-3 FTCA	6.11498	30	PFHxA				

HALO
PFAS Solutions

INNOVATION YOU CAN TRUST.
PERFORMANCE YOU CAN RELY ON.

HALO® and Fused Core® are registered trademarks of Advanced Materials Technology.

EPA 1633: 40 PFAS Calibration Standard




TEST CONDITIONS:

Column: HALO 90Å PFAS 2.7µm, 2.1 x 50 mm
Part Number: 92812-413
Delay Column: HALO® PFAS Delay, 2.7µm, 3.0 x 50 mm
Part Number: 92113-415
Mobile Phase A: 2 mM Ammonium Acetate in 95:5 water/acetonitrile (v/v)
Mobile Phase B: Acetonitrile
Gradient:


Time	%B	Flow Rate (mL/min)
0.0	2.0	0.35
0.2	2.0	0.35
4.0	30.0	0.40
7.0	55.0	0.40
9.0	75.0	0.40
10.0	95.0	0.40
10.4	2.0	0.40
11.8	2.0	0.40
12.0	2.0	0.35

Pressure: 195 bar
Temperature: 40 °C
Injection Volume: 2.0 µL
Sample Solvent: 96:4 Methanol/Water
MS System: Agilent 6495C
LC System: Agilent 1290Infinity II
Detector: MS (ESI-), MRM
Data Courtesy of: Lara Rosenberger



Analyte	Abbreviation	Retention Time (min)
Perfluorobutanoic acid	PFBA	1.79
Perfluoro-3-methoxypropanoic acid	PFMPA	2.51
3-perfluoropropyl propanoic acid	3:3 FTCA	3.57
Perfluoropentanoic acid	PFPeA	3.70
Perfluorobutanesulfonic acid	PFBS	4.12
Perfluoro-4-methoxybutanoic acid	PFMBA	4.21
Perfluoro(2-ethoxyethane)sulfonic acid	PFEESA	4.68
Nonafluoro-3,6-dioxahexanoic acid	NFDHA	4.94
1H, 1H, 2H, 2H-Perfluorohexane sulfonic acid	4:2FTS	5.03
Perfluorohexanoic acid	PFHxA	5.14
Perfluoropentanesulfonic acid	PFPeS	5.33
Hexafluoropropylene oxide dimer acid	HFPO-DA	5.45
Perfluoroheptanoic acid	PFHpA	6.10
Perfluorohexanesulfonic acid	PFHxS	6.18
2H,2H,3H,3H-Perfluorooctanoic acid	5:3 FTCA	6.18
4,8-dioxa-3H-perfluorononanoic acid	ADONA	6.21
1H, 1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2FTS	6.89
Perfluorooctanoic acid	PFOA	6.93
Perfluoroheptanesulfonic acid	PFHpS	6.99
Perfluorononanoic acid	PFNA	7.79
Perfluorooctanesulfonic acid	PFOS	7.81
3-perfluoroheptyl propanoic acid	7:3 FTCA	8.00
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid	9Cl-PF3ONS	8.28
1H, 1H, 2H, 2H-Perfluorodecane sulfonic acid	8:2FTS	8.59
Perfluorononanesulfonic acid	PFNS	8.61
Perfluorodecanoic acid	PFDA	8.61
N-methyl perfluorooctanesulfonamidoacetic acid	NMeFOSAA	8.98
Perfluorooctanesulfonamide	PFOSA	9.04
Perfluorodecanesulfonic acid	PFDS	9.35
Perfluoroundecanoic acid	PFUnA	9.37
N-ethyl perfluorooctanesulfonamidoacetic acid	NEtFOSAA	9.37
11-Chloroicosadecafluoro-3-oxaundecane-1-sulfonic acid	11Cl-PF3OUdS	9.74
Perfluorododecanoic acid	PFDoA	10.05
N-Methyl Perfluorooctanesulfonamide	NMeFOSA	10.51
N-Methyl Perfluorooctanesulfonamidoethanol	NMeFOSE	10.54
Perfluorododecanesulfonic acid	PFDoS	10.62
Perfluorotridecanoic acid	PFTriDA	10.66
N-Ethyl Perfluorooctanesulfonamidoethanol	NEtFOSE	11.00
N-Ethyl Perfluorooctanesulfonamide	NEtFOSA	11.00
Perfluorotetradecanoic acid	PFTA	11.19

Field Soil Sample



Analyte	Abbreviation	Retention Time (min)
Perfluorobutanoic acid	PFBA	1.79
Perfluoropentanoic acid	PFPeA	3.70
Perfluorobutanesulfonic acid	PFBS	4.12
Nonafluoro-3,6-dioxaheptanoic acid	NFDHA	4.94
Perfluorohexanoic acid	PFHxA	5.14
Perfluoroheptanoic acid	PFHpA	6.10
Perfluorohexanesulfonic acid	PFHxS	6.18
2H,2H,3H,3H-Perfluorooctanoic acid	5:3 FTCA	6.18
1H, 1H, 2H, 2H-Perfluorooctane sulfonic acid	6:2FTS	6.89
Perfluorononanoic acid	PFNA	7.79
Perfluorooctanesulfonic acid	PFOS	7.81
Perfluorodecanoic acid	PFDA	8.61
N-methyl perfluorooctanesulfonamidoacetic acid	NMeFOSAA	8.98
Perfluorooctanesulfonamide	PFOSA	9.04
Perfluoroundecanoic acid	PFUnA	9.37
N-ethyl perfluorooctanesulfonamidoacetic acid	NEtFOSAA	9.37
Perfluorododecanoic acid	PFDoA	10.05


Exemplary Ultra Short-Chain PFAS (USC)

Challenges:


- Retention
- Sensitivity
- Peak Shape

TFA

114.02 MW




C1



PFMtS

150.08 MW

C2




PFEtS


200.08 MW

PFPrA

164.03 MW



C3




PFPrS


250.09 MW

PFBA

214.04 MW



C4



PFBS

300.10 MW



Ultrashort-chain-PFAS (C1-C4)
DRE-A30000064MW

- Fluorine inductive effect causes low pKa's for compounds with an adjacent R-OH, R-COOH or R-SO₃H.
- Many compounds of interest are fully ionized under typical LC separation and ESI conditions.

Method Development USC PFAS


HILIC vs. Reverse Phase: Charged Surface Materials



HILIC: Charged Surface

EXAMPLE SEPARATION OF PFAS on HALO® HILIC-Plus Prototype Column

Conditions: 2.1 ID x 100 mm HALO HILIC-Plus Prototype, 0.5 mL/min; 35°C
85-65% B in 5 minutes, hold for 1 min., return to 85%;
B- 80% AcN/20%IPA/0.1% FA, A - 0.1% FA/10 mM AF in water




Rt and Compounds:

1.40 - 11-Chloroeicosafuoro-3-oxanonane-1-sulfonate	2.52 - Perfluorohexanoic acid
1.51 - 9-Chlorohexadecafluoro-3-oxanonane-1-sulfonate	2.76 - Perfluoropentanoic acid
1.56 - Perfluorooctanesulfonic acid	3.13 - Perfluorobutanoic acid
1.70 - Perfluorohexanesulfonic acid	3.76 - Perfluoropropionic acid
1.91 - Perfluorobutanesulfonic acid	5.38 - Trifluoroacetic acid
2.21 - Perfluorooctanoic acid	

Reverse Phase: Charged Surface

EXAMPLE SEPARATION OF PFAS on HALO® PCS-Phenyl-Hexyl

Conditions: 2.1 ID x 100 HALO PCS Phenyl Hexyl, 0.4 mL/min; 40°C
2-95% B in 7 minutes, hold for 1 min., return to 2%;
B- MeOH, A-5mM Amm. Formate, 0.05% FA in water
Thermo QE-HF in PRM mode






- | | |
|----------------------------------|----------------------------------|
| 1. Trifluoroacetic acid | 5. Perfluorobutanoic acid |
| 2. Perfluoromethanesulfonic acid | 6. Perfluoropropanesulfonic acid |
| 3. Perfluoropropionic acid | 7. Perfluorobutanesulfonic acid |
| 4. Perfluoroethanesulfonic acid | |

Method Development: Ultrashort Chain to Long Chain


Column Selection

- Reverse Phase vs. HILIC?
- Stationary Phase – Mixed
- Column Dimension
- SPP vs. FPP



Mobile Phase Optimization





- Acidic Modifier/ Concentration
- Buffer Modifier/ Concentration
- pH, Temperature



Dry Lab Optimization


Method Development: USC to LC

Stationary Phase, Mobile Phase, Temperature, Gradient




Aqueous Additive	Concentration
Formic Acid	0.01-0.1%
Acetic Acid	0.01-0.1%
Ammonium Formate	5mM
Ammonium Acetate	5mM
Ammonium Formate/ Formic Acid	5mM, 0.05%

Temperature
30°C
40°C
Organic Modifier
Acetonitrile
Methanol



Dry Lab®

HALO® PFAS PCS Phenyl-Hexyl: Separation of Ultrashort PFAS




Testing Conditions:

Mobile Phase: A: 5mM Ammonium Formate, 0.05% FA
B: Methanol
Gradient: 2-95%B in 7min
Instrument: Shimadzu 8060 NX
Injection: 2 µl (12.5µg/mL)
Temperature: 40 °C
Flow Rate: 0.4 mL/min.
Column: HALO 90Å PCS Phenyl Hexyl, 2.7µm, 2.1 x 100mm
Delay Column: HALO 160Å PFAS Delay, 2.7µm, 3.0x50mm




Gradient Conditions: 2-95% B in 7min, 2 min. hold, Flow Rate 0.4 mL/min




Aqueous Component Comparisons

PFAS Retention Time Shift in Formic Acid




- Increase USC retention with good peak shape including elution of LC PFAS
- Minimize interference with MS signal
- Common reagents in use or similar to previous methods
- Maintain full system pressure below 500 bar



- TFA improved distance from V_0
- Peak shapes are otherwise acceptable, although PFDoDeA somewhat broadened

PFAS PCS Phenyl-Hexyl un-normalized



1. Trifluoro acetic acid
2. Trifluoromethane sulfonic acid
3. pentafluoropropionic acid
4. pentafluoroethanesulfonic acid
5. heptafluorobutyric acid
6. perfluoropropanesulfonic acid

7. perfluoropentanoic acid
8. perfluorobutanesulfonic acid
9. perfluorohexanoic acid
10. perfluoroheptanoic acid
11. perfluorohexanesulfonic acid
12. perfluorooctanoic acid

13. perfluorononanoic acid
14. perfluorooctanesulfonic acid
15. perfluorodecanoic acid
16. perfluoroundecanoic acid
17. perfluorodecanesulfonic acid
18. perfluorododecanoic acid

19. perfluorotridecanoic acid
20. perfluorotetradecanoic acid
21. perfluorohexadecanoic acid
22. perfluorooctadecanoic acid

- NEW -



METHOD 534: DETERMINATION OF SELECTED PER- AND POLYFLUOROALKYL SUBSTANCES IN DRINKING WATER BY SOLVENT DILUTION AND LIQUID CHROMATOGRAPHY/TANDEM MASS SPECTROMETRY

1.1 Analyte List

Analyte ^a	Abbreviation	CASRN
Hexafluoropropylene oxide dimer acid ^b	HFPO-DA	13252-13-6 ^b
Perfluorobutanesulfonic acid	PFBS	375-73-5
Perfluorodecanoic acid	PFDA	335-76-2
Perfluoroheptanoic acid	PFHpA	375-85-9
Perfluorohexanesulfonic acid	PFHxS	355-46-4
Perfluorohexanoic acid	PFHxA	307-24-4
Perfluorononanoic acid	PFNA	375-95-1
Perfluorooctanesulfonic acid	PFOS	1763-23-1
Perfluorooctanoic acid	PFOA	335-67-1

2 Method Summary

A 3–6 mL drinking water sample is collected in a tared polypropylene tube. Isotopically substituted analogues are added as internal standards. The sample is dechlorinated in the laboratory by the addition 10 µg ascorbic acid. The sample is diluted with methanol to a final concentration of 50% methanol/50% aqueous, then 10 µL acetic acid is added. The sample is mixed, and an aliquot is transferred to an autosampler vial. The diluted sample is analyzed by liquid chromatography tandem mass spectrometry (LC-MS/MS) using a mixed-mode column with high retentivity for anions. Analytes and their internal standards are detected in the MRM mode. The concentration of each analyte is calculated using internal standard calibration.

17 Tables, Figures and Method Performance Data

Table 1. HPLC Method Conditions^a

Time (min)	% 2 mM aq. Ammonium acetate	% 2 mM ammonium acetate in methanol
Initial	100	0
2.0	100	0
3.0	50	50
15.0	5	95
17.0	5	95
18	100	0
20	100	0

^a Atlantis Premier BEH C18 AX VanGuard FIT Column, 2.5 µm, 2.1 x 50 mm, column temperature 50 °C. Flow rate of 0.4 mL/min; 200 µL injection into a 250 µL loop (while the loop was set to inject 200 µL, due to specific characteristics of the instrument sample loop the actual injected volume was 199 µL). The chromatogram in Figure 1 was obtained under these conditions.

Table 2. ESI-MS Method Conditions

ESI Conditions for Waters Absolute	
Polarity	Negative ion
Capillary needle voltage	0.90 kV
Cone gas flow	150 L/hour
Nitrogen desolvation gas	1000 L/hour
Desolvation gas temperature	400 °C

USC PFAS in Well Water ~10-100 ng/L



Testing Conditions:

Mobile Phase: A: 5 mM AF, 0.05% FA

B: MeOH

Gradient: 2-95%B in 7 min, 2 min hold



Injection: 20 μ L, 0.05% HOAc

Temperature: 40°C


Flow Rate: 0.5 mL/min

Column: HALO 90Å PCS Phenyl Hexyl, 2.7 μ m, 2.1 x 100 mm

Instrument: Agilent 1290 UPLC, 6495C QQQ




Short and Long PFAS in Soil Extract ~0.1-1 ng/g




Looking Forward: More!

- Capillary work/ Sensitivity Gains
 - Scale internal diameter (1.5, 1.0 mm, etc.)
 - Responding to an emerging Public Health/Regulatory need
 - Increased analysis demand will favor decreasing analysis time, solvent use, waste generation and overall cost
- Sample Prep
 - Isotopic Standards are great, but a real need is automating or eliminating sample preparation and quantification
 - Either large samples, or effective concentration, but Simpler and Cheaper
- Materials Development to enable higher sample throughput via mixed mode and multimodal Workflows:
 - One size fits all is appealing, but potentially not realistic. HILIC/AEX and HILIC/AEX with 2D approaches may offer the highest sensitivity, but a more complex workflow
- More Analytes, More Matrices: The more we look, the more we find!

ex. Bis(trifluoromethanesulfonyl)imide: used as an electrolyte in lithium-ion batteries and is showing up especially in landfill leachate



Staci L. Capozzi, Chunjie Xia, Matthew Shuwal, Gillian Zaharias Miller, Jeff Gearhart, Erica Bloom, Lennart Gehrenkemper, Marta Venier,
From watersheds to dinner plates: Evaluating PFAS exposure through fish consumption in Southeast Michigan, Chemosphere, Volume 345, 2023.




Conclusions/ Acknowledgements



- Advanced Materials Technology

- Superficially Porous Particle Technology
- Current EPA Methods in place and challenges for the Ultrashort Chain
- HALO PFAS PCS Phenyl Hexyl



- Special Thanks to...



- Technical Resources/ Support

- halocolumns.com







Image credit: NIEHS.

Questions?



Technical and Marketing Materials:

- www.halocolumns.com

Technical Support:

- support@advanced-materials-tech.com

A screenshot of the "TECHNICAL RESOURCES" website. The header is dark blue with the text "TECHNICAL RESOURCES" in white. Below the header, there are five circular icons representing different resource types: Application Notes (document with pencil), Conference Papers (person at a presentation board), Product Literature (open book with magnifying glass), Technical Documents (stack of papers), and Videos (computer monitor with play button). Each icon has its name written below it. Below these icons is a search and filter section. It includes a text input field labeled "Enter Keyword" with the value "PFAS". To the right of this are four dropdown menus: "Market Segment", "Particle Size", "Pore Size", "Product Segment", "Separation Mode", and "USP Number". At the bottom of this section are two buttons: "CLEAR FILTERS" and "APPLY FILTERS".

Conner McHale
Technical Support Specialist
Advanced Materials Technology

cmchale@advanced-materials-tech.com
Phone: 1-302-992-8060 *1124



advancedmaterialtechnology



halocolumns.com



3521 Silverside Road, Suite 1-K
Quillen Building
Wilmington, DE 19810



(302) 992-8060

